

Envirostar Sharps Containers

For Non-Infectious, Non-Biohazard Sharps

Keeping noninfectious sharps separate increases safety and cuts cost. The EnviroStar™ non-infectious waste sharps disposal container, help handle broken light bulbs, slides, glassware, and other sharp materials that can cause cuts but not blood borne infection. The containers are manufactured with 100% recycled polypropylene and are uniquely colored to easily distinguish proper waste segregation of non-infectious sharps. The EnviroStar™ containers are autoclavable and incinerable, emitting only carbon dioxide and water vapors when burned. Containers meet or exceed all JCAHO and OSHA requirements.



American MasterTech
scientific laboratory supplies

Size:

2 Gallon
8 Gallon
18 Gallon

Dimensions:

10"H x 7.25"D x 10.5"W
17.75"H x 11"W x 15.5"D
26"H x 12.75"D x 18.25"W

Item#:

CPW02501E
CPW02503E
CPW02505E

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.