



American MasterTech
scientific laboratory supplies

Fan Pad GL™

NASA's #1 Choice for formaldehyde neutralizing pads!

FAN PAD GL™ super absorbent pads are specially treated with Formalex® for the ultimate destruction of formaldehyde, 10% formalin, and glutaraldehyde solutions. FAN PAD GL™ also eliminates the hazardous vapors associated with these poisonous solutions.

FAN PAD GL™ are perfect for use in: Histology, Pathology, Research, Small Spills, Tissue Transport, Grossing Examination, Operating Room, Formalin storage, Tissue storage, Gross Dissection, Labor/Delivery, Embalming suites, Dermatology Labs, Autopsy & Morgues, G.I. Lab & Endoscopy, and Biology & Anatomy.

Mini Fan Pad absorbs 20-30ml of 10% Formalin

Size: 8" x 11"

Qty: Case(330)

Item#: CWFAN-55

Small Fan Pad absorbs and neutralizes up to 300cc of 10% Formalin

Size: 6" x 14"

Qty: Case(48)

Item#: CWFAN-614

Large Fan Pad absorbs and neutralizes up to 900cc of 10% Formalin

Size: 8" x 25"

Qty: Case(30)

Item#: CWFAN-825



Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.