

# POLY/Fin™ Paraffin . . .

*Does Your Paraffin measure up?*

- ◆ *Minimal compression & cell disruption*
- ◆ *Can safely be used in P/V processors*
- ◆ *Ideal handling characteristics*
- ◆ *Cut < 2 microns without icing*



**POLY/Fin™** is a mixture of highly refined paraffins and co-polymers in a convenient wax pellet. This unique formulation provides optimal tissue support while maintaining exceptional clarity. TBS' **POLY/Fin™** has been institutionally tested and found to function flawlessly in pressure/vacuum fluid flow processors. Low (55° C) melting point eliminates any tissue distortion caused by excessive heat during processing.

American MasterTech  
scientific laboratory supplies

1330 Thurman Street  
Lodi, CA 95240

Phone: 800.860.4073  
Fax: 209.368.4136

[www.americanmastertech.com](http://www.americanmastertech.com)