

Histology Innovation for a NEW Generation

Embedding

Forceps Warmer



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

- ✔ **Three Convenient Wells Provides Space for Multiple Forceps**
The three wells allow for multiple forceps to be used thus eliminating the need to wait on individual forceps heating, thus increasing lab efficiency.
- ✔ **Consistent Temperature Prevents Paraffin Wax Build-Up**
Placing the forceps into one of the three wells heated to 80°C the paraffin wax is melted off and allows easier tissue orientation at embedding centers.
- ✔ **Angled Forceps Wells Makes Retrieving & Replacing the Forceps Easy**
- ✔ **Portable for Easy Addition to Embedding Station**
The electric Forceps Warmer can be placed on top of any embedding station due to its small compact size.
- ✔ **Switch Protector Casing Prevents Paraffin Wax Contamination**
The 2 piece design is durable, safe and resistant to paraffin contamination.