

TEC-III™

Tissue Embedding Center



■ Dispensing Console Creates Versatility

An adjustable dispenser paddle accommodates mega cassettes and provides greater workspace, and two built-in paraffin block trimmers eliminate the need for extra equipment. White LED illumination creates non-glare visualization for superior embedding. Eight heated forceps well ports keep the workspace clean.

■ Large Capacity Thermal Console

The thermal console will hold baskets from any processor, and has capacity for two baskets (200 cassettes) or three baskets (300 cassettes) when trays are removed. The high profile cover opens and closes easily to provide quick access to cassette storage.

■ Cryo Console Work Surface Accommodates More Cassettes/Base Molds

The cryo console has a large cooling surface that holds up to 60 cassettes/molds. Variable temperature settings can be set between 0 and -12 degrees C. Steady temperature distribution provides consistency over the entire console surface. This unit can be placed to the right or the left of the dispensing console for workflow flexibility.

Protect every life story™ in your lab.

Product Specifications:

| DIMENSIONS: | | DISPENSING CONSOLE ADJUSTABLE TEMPERATURE RANGE: | |
|---------------------------------|-------------------|--|--------------|
| Dispensing Console (W x D x H): | 13 x 25 x 16 (in) | Paraffin Reservoir: | 40°C to 70°C |
| Thermal Console (W x D x H): | 13 x 24 x 15 (in) | Dispensing Nozzle: | 40°C to 70°C |
| Cryo Console (W x D x H): | 13 x 24 x 15 (in) | Thermal Chambers: | 40°C to 70°C |
| Total Unit (W x D x H): | 39 x 25 x 16 (in) | Auxiliary Heated Forceps: | 60°C to 75°C |
| CAPACITY: | | CRYO CONSOLE ADJUSTABLE TEMPERATURE RANGE: | |
| Paraffin Reservoir: | 5 Liters | Adjustable Range: | 0°C to -12°C |

| Catalog # | Description |
|-------------|--|
| TEC 2130-CU | TEC-III Embedding Center, includes: <ul style="list-style-type: none">• Dispensing console• Thermal console• Cryo console• Power cord• User manual |

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2017 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com