



FISHER ISOTEMP® STANDARD OVENS

IDEAL FOR GENERAL LABORATORY DRYING

ISTOTEMP® STANDARD OVEN DESCRIPTION:

- Temperature range 50 to 225°C (122 to 437°F)
- Three sizes available
- Gravity-flow convection
- Electronic control
- Highest quality economical oven available
- Quick and easy temperature setting
- Fast temperature recovery time
- Adjustable shelf heights
- Exterior remains at safe temperature
- Back-up temperature control
- Indicator lights show which thermostat is controlling oven temperature
- Durable powder-coated steel exterior
- Corrosion resistant, easy to clean, stainless-steel interior
- Oven temperature can be monitored in port on top
- Doors open 180°
- Door features two grabber-type latches
- Non-slip black synthetic rubber feet • 120V models are CSA approved

ISTOTEMP® STANDARD OVEN SPECIFICATIONS:

Item#:	Chamber Dimensions Inches (cm)			Exterior Dimensions Inches (cm)			Volume Cubic ft. (L)	Air Changes Per Hr.	Recovery Time	Control Sensitivity	Shipping Weight Lb. (kg)
	D	W	H	D	W	H					
EQO0102	18 (46)	18 (46)	13½ (34)	23½ (60)	25.8 (65)	26 (66)	2½ (70.8)	24	2 min.	±0.5°C	115 (52)
EQO0102	18 (46)	18 (46)	20 (51)	23½ (60)	25.8 (65)	33 (84)	3¾(106.2)	16	4 min.	±0.5°C	125 (57)
EQO0103	18 (46)	18 (46)	26½ (67)	23½ (60)	25.8 (65)	39.3 (100)	5.0(141.6)	12	4 min.	±0.5°C	140 (64)

Copyright ©2012-2018 American MasterTech Scientific Incorporated. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.