

Histology Innovation for a NEW Generation

Sectioning/Incubation

SHURDry™ Slide Dryer III



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

- ✔ **See-through Glass Allows for Easy Content Monitoring**
Is convenient for cleaning and provides the histologist the visibility needed for examining contents without releasing heat by opening doors.
- ✔ **Keeps Paraffin Liquefied in Between Runs**
SHURDry™ III design makes storage of melted paraffin for microwave processing highly efficient. It eliminates the need to remelt the paraffin and can save on waste & consumable cost.
- ✔ **Large Capacity Interior Accommodates An Assortment Glassware & Racks**
Adjustable shelf can accommodate any size slide rack or bowl that needs to be kept warm.
- ✔ **Operator Defined Programs to Control Drying Times**
The control panel utilizes six function keys that are available with two built-in timers for multi-user drying times. Keys can be programmed to a different time and an alarm will sound when an item is ready to be removed.
- ✔ **Small Footprint Fits on Counter Tops**