

Paraffin Dispenser PRO™

Paraffin Dispenser



- **Large Reservoir Keeps Liquid Paraffin Ready**
Liquid paraffin will be ready when you need it with the 6.25 gallon reservoir.
- **Controls Keep Paraffin at a Constant Temperature**
Two heating elements work with a microprocessor to keep paraffin within 1°C of desired temperature. When new paraffin pellets are added to the reservoir, sensors detect the change in temperature and automatically adjust. Overheating protection assures safety.
- **Accommodates a Variety of Large Bowls/Dishes**
High spigot and concave cut-out allow users to conveniently place and maneuver a variety of large paraffin collection bowls / dishes for added flexibility.
- **Easy Paraffin Viewing Without Contamination**
The see-through Lexan lid offers effortless level monitoring. This unique access allows the user to refill paraffin pellets without lid removal.

Protect every life story™ in your lab.

Product Specifications:

Dimensions: WxDxH (in) [cm] (14x15x23) [35.6x38.1x58.4]

Weight: (lbs) [kg] (45) [20.5]

Voltage: 100-120VAC, 50/60 Hz, 4.5A and 220-240 VAC 50/60Hz, 2.4A

Certification: TUV

Catalog #	Description
PDP-120	Paraffin Dispenser, microprocessor controlled, 6.25 gallon capacity w/heated non-drip spigot, lexan lid, built-in clock, ambient to 65°C. UL/CSA/CE; Specifications: 100-120VAC, 50/60Hz, 4.5A WxDxH: (14x15x23)in, (35.6x38.1x58.4)cm; Weight: 45lbs, 20.5kg
PDP-220	Paraffin Dispenser, microprocessor controlled, 6.25 gallon capacity w/heated non-drip spigot, lexan lid, built-in clock, ambient to 65°C. UL/CSA/CE; Specifications: 220-240VAC, 50/60Hz, 2.4A; WxDxH: (14x15x23)in, (35.6x38.1x58.4)cm; Weight: 45lbs, 20.5kg

More Information

844.643.1129

www.general-data.com/hc

General Data Healthcare helps labs enhance productivity and improve workflow with innovative products, solutions and support that focus on delivering better patient safety & care.

©2018 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
4043 McMann Road, Cincinnati, OH 45245
844.643.1129

Email: talktous@general-data.com
Web: www.general-data.com/healthcare

Service:
PH: 800.245.5029
Email: customersupport@general-data.com