

Histology Innovation for a NEW Generation

Sectioning/Incubation

SHURDry™ Slide Dryer II



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

- ✓ **See-through Glass Allows for Easy Monitoring**
Convenient for cleaning and provides the histologist the visibility needed for examining contents without releasing heat by opening doors.
- ✓ **Perfect Solution for Small Volume Labs**
With an internal dimension of 10 x 4.75 x 3.25 (in) it can accommodate multiple slide stainer racks.
- ✓ **Easy-To-Read Indication Lights Assist With Heat Control**
The Amber indication light located by the temperature control, informs the user when the heat is on. The Red indicator lights offers protection, and informs the user if overheating occurs.
- ✓ **Dries Slides Quickly**
The Slide Dryer II is able to dry two slide stainer racks within 20 minutes to 1 hour.
- ✓ **Small Footprint Fits on Counter Tops**

Dimensions: WxDxH (in)(cm)	(11 x 10 x 5) (26 x 28 x 13) (in) (cm)
Weight :lbs(kg)	11 (5)
Voltage:	100-120 VAC, 50/60Hz, 2A and 220-240VAC, 50/60Hz, 1A


American MasterTech
scientific laboratory supplies

1330 Thurman Street
Lodi, CA 95240

Phone: 800.860.4073
Fax: 209.368.4136

www.americanmastertech.com

©2013 General Data Healthcare, Inc. All rights reserved. Specifications are current at the time of publication; however, they are subject to change without notification. All trademarks and registered trademarks are property of General Data Healthcare, Inc., Cincinnati, Ohio, USA.



General Data Healthcare, Inc.
Headquarters: 4354 Ferguson Drive, Cincinnati, OH 45245
TBS Facility: 3014 Croasdaile Drive, Durham, NC 27705

Sales Inquiries:
PH: 844.643.1129
Email: talktous@general-data.com
Web: www.general-data.com/hc

Service:
PH: 800.245.5029
Email: customersupport@general-data.com