

Histology Innovation for a NEW Generation

Sectioning/Microtomy

HISTO**O**rientator™ Specimen De-Wrinkle



**Proven Technology. Superior Reliability.
Protecting Every Life Story In Your Lab.**

- ✓ **De-Wrinkles & Removes Bubbles on Paraffin Sections**
The constant warming temperature at 70°C will quickly and easily heat slides to the necessary temperatures to remove wrinkles and bubbles on your sections.
- ✓ **Rapid Attachment of Fresh Frozen Sections**
HISTO**O**rientator™ can be kept on your cryostat to ensure solid attachment of frozen sections to slides to prevent loosening of the tissue during the staining agitation.
- ✓ **Requires Limited Laboratory Space**
Compact size will not take up valuable laboratory space. The 2"x 2" (5.1x5.1cm) pad allows two slides to be placed on the plate at the same time.
- ✓ **Light Indicator for ON/OFF Safety**
LED Select Switch indicates whether the HISTO**O**rientator™ is the ON/OFF mode.