

SlideFile™ Storage System



- Holds up to 400 slides
- Stackable for space efficient storage
- Available in 5 distinctive colors for easy color coding
- Read bar codes without having to remove slides from the box



This impact resistant SlideFile™ storage system can hold up to 400 slides (3 x1in) per unit and is stackable for space efficient storage. The main feature to the SlideFile™ system is a removable tray inside the storage box having 100 individual numbered slots. For space saving purposes, you can double the amount of slides simply by storing 2 slides per slot. And for maximum storage space, simply remove the tray and line up 400 slides in 3 rows. Not autoclavable. Dimensions: 82 mm x 245 mm x 86 mm.



American MasterTech
scientific laboratory supplies

Color:	Item#:
White	FMS00640E
Blue	FMS00604E
Green	FMS00613E
Pink	FMS00632E
Yellow	FMS00641E

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.