


HELDSECURE™

Slide & Cassette Block Storage

A 10% savings in comparison to similar storage solutions!



Store 2,000 slides or 500 cassettes per box!

Looking for an economical way to archive your cassettes and slides? The HeldSecure Slide/Cassette Storage Box is a durable, stackable, organizational system for archiving your slides and cassettes. They are made from quality, corrugated, non-coated fiber board and come in a case of 10 unassembled boxes. The HeldSecure Slide Storage Box includes 50 In/Out cards, 25 Index Cards, and four double row trays that can accommodate up to 2,000 slides. The HeldSecure Cassette Block Storage Box includes 50 In/Out cards and eight double row trays that can accommodate up to 500 plus cassette blocks. Additional cards and tray labels available. Both the cassette and slide boxes are stackable with the H400 & H800 Tissue Block/Slide filing system.



American MasterTech
scientific laboratory supplies

Description:

HeldSecure Slide Storage / Case (10)
HeldSecure Cassette Storage / Case (10)
IN-OUT Cards for Slide Box / Pack (1000)
IN-OUT Cards Cassette Box/ Pack (1000)
Index Cards for Slide Box / Pack (500)
Removable Tray Labels / Roll (500)

Item#:

FMS1401C
FMS1402C
FMS1403P
FMS1405P
FMS1404P
FMS1406P

Copyright © 2012-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.