



Looking for a convenient way to properly collect, transport, and identify a prostate biopsy?

PROSTATE Biopsy Kits

- › Pre-Filled with 10% Formalin
- › Pre-labeled, positioned vials on a location diagram
- › Additional, pre-printed labels
- › Specimen Transport Bag Included

Prostate Biopsy Kits are a convenient way to properly collect, transport, and identify a prostate biopsy. These Prostate Biopsy Kits are very simple to use. Just collect the appropriate tissue and place it into the pre-labeled, positioned vials located on a corrugated insert. The corrugated insert is designed in a way so health care professionals can easily identify where the tissue was taken from. A set of additional labels allows duplicate identification of the specimen biopsy location so that when the vials are removed from the insert, the location of the tissue from the prostate can still be easily identified.

The Prostate Biopsy Kits come in six, eight, or twelve 20ml vials that are pre-filled with 10ml of 10% formalin. Each vial has been pressure-tested according to the 95kPa standards for leak resistant transport by air or ground. The Prostate Biopsy Kit also includes absorbent material and a specimen transport bag to secure the vials and corrugated insert during transport.

The 6 vial Prostate Biopsy Kits come in a case of 24 kits, the 8/12 vials Prostate Biopsy Kits come in a case of 30 kits, and the 16 vial Prostate Biopsy Kits come in a case of 20 kits.

Description:	Qty:	Item#:
6 Vial Prostate Biopsy Kit	24 Kits	FXF0206C
8 Vial Prostate Biopsy Kit	30 Kits	FXF0208C
12 Vial Prostate Biopsy Kit	30 Kits	FXF0212C
16 Vial Prostate Biopsy Kit	20 Kits	FXF0216C

