



ColorFrost Slides



American MasterTech
scientific laboratory supplies

ColorFrost™ Slides are super clean and have a smooth colored writing surface that is designed not to smear important information. The smooth colored tabs make identifying and categorizing slides easy and the tabs are resistant to common laboratory, chemicals and reagents. The positive charge on the adhesive slide allows for maximum adhesion, which eliminates time consuming coating steps in the laboratory. Plus slides also virtually eliminate background staining in standard H&E stains. 10 gross in a case.

COLOR:	Regular / Non-Adhesive Slides (10 gross) Item#:	Plus/ Adhesive Slides (10 gross) Item#:
Aqua	GL9951A	GL9951PA
Blue	GL5951	GL5951PB
Burnt Orange	GL9951BO	GL9951PBO
Gold	GL9951GL	GL9951PGL
Green	GL6951	GL6951PG
Lilac	GL9951L	GL9951PL
Pink	GL7951	GL7951PP
Tan	GL9951T	GL9951PT
Yellow	GL8951	GL8951PY

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.