



American MasterTech
scientific laboratory supplies

POST OFFICE BOX 2539 LODI, CALIFORNIA 95241 TELEPHONE: 1 (800) 860 4073 FACSIMILE: 1 (209) 368 4136

Adenna Gold Latex

Powder Free Examination Gloves



Adenna Gold Latex Powder Free Examination Gloves are the ultimate choice of dental and medical professionals who demand uncompromising quality and protection, with the minimum risk of potential glove related chemical or protein allergies. Specially manufactured with light chlorination rinse, these gloves maintain a non-slip, 'tacky' surface that provides excellent grip for both dry and wet applications.

Features & Benefits

- Powder Free minimizes the risk of developing latex protein or chemical residue related allergies;
- Eliminate powder related contamination in critical exam or in clean room environment;
- Full textured surface provides superior grip in wet or dry applications;
- Smooth interior finish for easy donning and comfortable wear.