



Monoclonal Antibodies to CD43

Technical Data Sheet

Reagent Category

Monoclonal primary antibody to CD43

Antibody Formats Supplied

0.2 ml of anti CD43, Concentrate

PRODUCT# IMI15852E (Concentrate)

SOURCE: Cultured supernatant

IMMUNOGLOBULIN ISOTYPE: IgG1

CLONE: MT1

REACTIVITY

CD43 reacts with myelocytic cells, it reacts with both cytoplasmic and membrane antigens that are expressed on myelocytic and T cells. This antibody is routinely used in identification of T cell lymphomas and normal T cells. This antibody is suitable for use on routinely fixed paraffin embedded and frozen sections. This antibody is compatible with formalin fixative as well as most commonly used histological fixatives.

POSITIVE CONTROL TISSUE: Mantle cell lymphoma, Tonsil

STAINING PATTERN: Membrane and cytoplasmic staining

Antibody FORMAT, DILUTION and INCUBATION

This affinity-purified monoclonal antibody is provided in phosphate buffered saline pH 7.2 containing bovine serum albumin (BSA) and 0.05% sodium azide.

0.2 ml of Concentrate, for paraffin section staining: Dilute 1:20 for 30 minute incubation. The dilution factor for frozen section staining and other applications are considerably higher and must be determined by the end user for the application and the secondary staining reagents employed. This antibody best stains when formalin fixed paraffin embedded sections are retrieved.

For Background-Free and easy retrieval, Innovex Uni-Trieve (60°-70°C) universal retrieval solution is highly recommended. To retrieve with Uni-Trieve: Warm up water bath to 60°-70°C and fill a slide holding container with Uni-Trieve solution. Place slides in a slide rack and immerse slides in a vessel or a coplin jar containing Uni-Trieve solution, bring the Uni-Trieve solution and slides to 60°-70°C and let slides incubate at 60°-70°C for 30 minutes. Rinse slides with 2 quick changes of water; No cooling period is required. Proceed with immunostaining.

THE SUGGESTED STARTING DILUTION FACTORS ARE MERE GUIDELINES ONLY. OPTIMAL DILUTION FACTORS MUST BE DETERMINED BY THE END USER FOR THE SECONDARY STAINING REAGENTS, SPECIMEN, STAINING PROTOCOL AND APPLICATION EMPLOYED.

CONTINUED OVER



American MasterTech

scientific laboratory supplies

POST OFFICE BOX 2539 LODI, CALIFORNIA 95241 TELEPHONE: 1 (800) 860 4073 FACSIMILE: 1 (209) 368 4136

The above dilutions and incubation times were derived by the use of Innovex STAT-Q and HISTO-STAT staining systems. For less sensitive detection systems these dilution factors may vary and they should be determined by the end user for specific detection system employed.

APPLICATIONS

This antibody is intended for use in immunolocalization of CD43 by a variety of immunoassays such as Immunochemistry, fluorescent labeling, ELISA, Immunoblotting and dot blotting. In immunochemistry, this antibody can be used in staining routinely fixed paraffin embedded or frozen histological sections for the purpose of qualitative localization of CD43. This antibody can also be used for quantitative flow cytometric assays in the indirect method.

SPECIES REACTIVITY

This antibody reactivity with human has been substantiated; other species reactivity may be characterized by the end user.

STORAGE CONDITIONS

Store in refrigerator at 2-8°C through expiration date noted on the vial.

IMPORTANT NOTE:

The interpretation of test results is the sole responsibility of the end user.

FOR PROFESSIONAL AND RESEARCH USE ONLY