



Polyclonal Antibody to Hepatitis B surface antigen (HBsAG)

Technical Data Sheet

Reagent Category

Polyclonal primary antibody to
Hepatitis B surface antigen

Antibody Formats Supplied

□ 0.2 ml of Hepatitis B surface Antigen (HBsAG), Concentrate

PRODUCT# IMI24652E (Concentrate)

SOURCE: Rabbit purified IgG

IMMUNOGEN: Hepatitis B envelope protein

POSITIVE CONTROL TISSUE: Infected liver

REACTIVITY

This antibody recognizes hepatitis B surface antigen (HBsAg) and it reacts with the outer envelope of surface protein of hepatitis B virus infective virion. This antibody is non-cross reactive with the core component of the infective domain for hepatitis B virus. This antibody is suitable for use on routinely fixed paraffin embedded and frozen sections. This antibody is compatible with formalin fixative as well as most commonly used histological fixatives.

Antibody FORMATS, DILUTIONS and INCUBATIONS

This affinity purified polyclonal anti human Hepatitis B surface antibody is provided in phosphate buffered saline pH 7.4 containing bovine serum albumin (BSA) and 0.1% sodium azide.

□ **0.2 ml of Concentrate; For 30-minute incubation, dilute 1:50-1:100.** The optimal incubation time the working titer should be determined by the end user. **This antibody best stains when formalin fixed paraffin embedded tissue sections are retrieved.**

For Background-Free and easy retrieval, Innovex Uni-Trieve (60°C-70°C) universal retrieval solution is highly recommended. To retrieve with Uni-Trieve: Warm up a water bath to 60°C-70°C and fill a slide holding container with Uni-Trieve solution. Place slides in a slide rack and immerse slides in a vessel or a coplin jar containing Uni-Trieve solution, bring the Uni-Trieve solution and slides to 60°C-70°C and let slides incubate at 60°C-70°C for 30 minutes. Rinse slides with 2 quick changes of water; No cooling period is required. Proceed with immunostaining.

The above dilutions and incubation times were derived by the use of Innovex STAT-Q and HISTO-STAT detection systems and the use of Enhancing wash buffer (see special usage guidance below). For less sensitive detection systems these dilution factors may vary and they should be determined by the end user for specific detection system employed.

APPLICATIONS

This antibody is intended for use in immunolocalization of Hepatitis B surface antigen (HBsAG) by a variety of immunoassays such as Immunocytochemistry, fluorescent labeling, ELISA, Immunoblotting and dot blotting. In immunocytochemistry, this antibody can be used in staining routinely fixed paraffin embedded or frozen histological sections of human tissues for the purpose of quantitative localization of hepatitis B surface antigen (HBsAG). This antibody can also be used for quantitative flow cytometric assays in the indirect method.

CONTINUED OVER



American MasterTech
scientific laboratory supplies

POST OFFICE BOX 2539 LODI, CALIFORNIA 95241 TELEPHONE: 1 (800) 860 4073 FACSIMILE: 1 (209) 368 4136

SPECIES REACTIVITY

This antibody reactivity with human tissues has been substantiated, other species reactivity may be characterized by the end user.

STORAGE CONDITIONS

Store in refrigerator at 2-8°C through expiration date noted on the vial.

IMPORTANT NOTE:

The interpretation of test results is the sole responsibility of the end user.

FOR PROFESSIONAL AND RESEARCH USE ONLY