


Capture dirt before it enters the workplace



ADHESIVE STICKY MATS



STICKY MATS 18"x 36"(240Sheets)

Item#: LWF0301

STICKY MATS 25"X 45"(240Sheets)

Item#: LWF0302

Adhesive sticky mats are designed to remove dirt, wax, and contaminants from your shoes when entering or leaving a work area. They are made with an applied antimicrobial agent to assure long lasting protection from mildew, bacteria, and mold. When a mat becomes dirty, just peel off the used mat sheet and a new sheet is ready for use. Numbered corner tabs ensure one-at-a-time sheet removal. Available in two sizes: 18" x 36" and 25" x 45".



American MasterTech
scientific laboratory supplies

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.