

STAINTRAY™



American MasterTech
scientific laboratory supplies

A User Friendly Approach to Immunohistochemistry by Simport.

STAINTRAY™ is not only suitable for routine staining requiring a humid chamber but is ideal for Hematology, Cytology and Microbiology laboratories. Manipulation is made safe and easy by using only one hand.



The STAINTRAY™ has a black base made of tough ABS plastic, which is resistant to a wide range of chemicals. The bigger STAINTRAY™ will accept up to 20 slides on 4 plastic rails and the smaller STAINTRAY™ will hold up to 10 slides on 2 plastic rails. The plastic rails are covered with a polymer strip for holding slides in place, even if tray is held at an angle.

* Avoid chlorinated hydrocarbons.

When humidity is needed, wells between the raised rails will hold water securely without splashing, not only to avoid water touching the slides but to make them more easily retrieveable. Excess staining solution dripping from the slides, drains into the base where it can be emptied with a drain plug. Four rubber feet ensure greater base stability. Units are stackable for space saving purposes.



Two covers are available:

- A clear one allowing for visual examination. Made of PETG with a temperature range of -20°C to 60°C.
- A black lid for fluorescent work. Made of ABS with a temperature range of -80°C to 80°C.

Dimensions:

STAINTRAY™ 20 Slide: 15" x 9³/₈" x 1³/₄"

STAINTRAY™ 10 Slide: 9³/₈" x 9³/₈" x 1³/₄".

DESCRIPTION:

StainTray 20 Slide with Clear Lid LWS11CL
 StainTray 20 Slide with Black Lid LWS11BK
 StainTray 10 Slide with Clear Lid LWS10CL
 StainTray 10 Slide with Black Lid LWS10BK

ITEM#:

DESCRIPTION:

StainTray 10 Slide Black Lid LWS1001BK
 StainTray 10 Slide Clear Lid LWS1001CL
 StainTray 20 Slide Black Lid LWS1101BK
 StainTray 20 Slide Clear Lid LWS1101CL

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.