

SCIOLOGEX MicroPette Pippettors

Accurately sample and dispense liquid volumes with the single, eight, or twelve channel pippettors.

The Micropette is an adjustable volume, easy to calibrate pipettor, from SCIOLOGEX. These pipettors are ergonomically designed and lightweight which helps reduce physical stress and discomfort during use. The tip ejector allows for easy one-handed operation, while the ejector collar and tip cone can be removed for easy cleaning and maintenance.

A Tip ejector allows convenient one-handed operation.

B Simply turn the plunger button for volume selection.


C Finger support keeps the instrument in place with minimum user effort.

D Clear and easy-to-read display with big numbers.

E Ejector collar and tip cone can be removed for easy cleaning and maintenance.


F Dispensing head rotates for optimal pipetting convenience.

G Durable tip cone materials provide excellent chemical resistance.



Single Channel Pippettor

Volume(μl):	Increments(μl):	Item#:
0.1-2.5	0.05	PPP0101
0.5-10	0.1	PPP0102
2-20	0.5	PPP0103
5-50	0.5	PPP0104
10-100	1	PPP0105
20-200	1	PPP0106
100-1000	1	PPP0107
1000-5000	5	PPP0108



Eight Channel Pippettor

Volume(μl):	Increments(μl):	Item#:
0.5-10	0.1	PPP0201
5-50	0.5	PPP0202
50-300	5	PPP0203

Twelve Channel Pippettor

Volume(μl):	Increments(μl):	Item#:
0.5-10	0.1	PPP0301
5-50	0.5	PPP0302
50-300	5	PPP0303



American MasterTech
scientific laboratory supplies

Copyright ©2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.