

CANCER

DIAGNOSTICS, INC

Margin Marking Dyes



American MasterTech
scientific laboratory supplies

2oz. Bottle

Highest quality margin marking dyes for MOHs techniques, serial sections, embedding orientation, etc. Each individual color will remain bright and permanently adhered to surgical specimens through water, alcohols, clearing agents, and fixatives. 100% Satisfaction Guaranteed!

Color:	Item#:
BLACK	STD00103E
BLUE	STD00104E
GREEN	STD00113E
ORANGE	STD00130E
RED	STD00134E
VIOLET	STD00139E
YELLOW	STD00141E

8oz. Bottle

Color:	Item#:
BLACK	STD00203E
BLUE	STD00204E
GREEN	STD00213E
ORANGE	STD00230E
RED	STD00234E
VIOLET	STD00239E
YELLOW	STD00241E

Squeezable Marking Dyes



Each box contains 15 single-colored, re-closable, squeezable dye applicators. Each applicator holds 4ml.

Color:	Item#:	Color:	Item#:
Black	STD00303P	Red	STD00334P
Blue	STD00304P	Violet	STD00339P
Green	STD00313P	Yellow	STD00341P
Orange	STD00330P		

Prefilled Squeezable Dye Kit

Kit consists of 18 multi-use, multi-colored, reclosable squeeze applicator marking dyes. Each applicator holds 4ml. Dye Kit contains 3 of each color.

Item#: STD00459P

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.