

Davidson®

Margin Marking Dyes



American MasterTech
scientific laboratory supplies



The Davidson Marking System® has been used for years by physicians around the world as an unequivocally reliable method for color identification of surgical margins. Developed by Dr. Terence M. Davidson, M.D., the Davidson® Marking Dyes are specifically designed to withstand immersion into high concentrations of alcohols, fixatives, and clearing agents.

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.

The Davidson Marking System®



Available in seven easy-to-see colors, the dyes are applied directly to excised tissue for the purpose of orientation, multiple biopsy identification, or surgical margins. The dyes can be applied to frozen sections, fixed tissue, and freshly excised specimens.

2oz

Color:	Item#:
Black	STD01003E
Blue	STD01004E
Green	STD01013E
Orange	STD01030E
Red	STD01034E
Violet	STD01039E
Yellow	STD01041E

8oz

Color:	Item#:
Black	STD01103E
Blue	STD01104E
Green	STD01113E
Orange	STD01130E
Red	STD01134E
Violet	STD01139E
Yellow	STD01141E

2oz KIT

Kit includes all 2oz. dyes except orange and violet.

Item#: STD01308P

Brush Caps



For use with Davidson® 2oz. Margin Marking Dyes. 10 in a pack.

Item#: STD01400P

3cc Mini Tissue Dye



Contains 40 ea. Davidson® 3cc Mini tissue dyes in your choice of green, yellow, black, red, blue, orange or violet. These dyes are small enough to be single-use or can be re-capped for multiple use

3ml. Pk(40)

Color:	Item#:	Color:	Item#:
Black	STD00803P	Red	STD00834P
Blue	STD00804P	Violet	STD00839P
Green	STD00813P	Yellow	STD00841P
Orange	STD00830P		

DyeRack

Universal Margin Marking Dye Rack



Isn't it frustrating buying a marking dye rack that only holds one size/brand? With the Universal DyeRack from American MasterTech you can now store different sizes and brands of marking dyes in one convenient place! These Dye Racks are constructed of white acrylic and contain eight built-in holes for marking dyes and two holes for wooden applicator sticks. The Universal Dye Rack holds both Davidson's and Cancer Diagnostic's 2 oz. and 8 oz. bottles, as well as MasterTech's 1 oz. and 4 oz. bottles.

Item#: STD10