

# MasterTech Margin Marking Dyes



American MasterTech  
scientific laboratory supplies

Each MasterTech Margin Marking Dye is designed to withstand Fixatives, Alcohols, and Clearing Agents. MasterTech Margin Marking Dyes are not stains, but special "paints", designed to coat the area where applied without staining cellular material which could be of diagnostic value.

## Now Available In Bright, Contrasting White!

Identifying a tissue margin is easy as ever with these new MasterTech White Marking Dyes. The white dye gives users enhanced visibility during embedding, block orientation, and tissue flotation. Margin appears bright white to the naked eye for fast orientation and dark gray to black under the microscope.



### 1/4oz

Color:	Item#:
Black	STD02003E
Blue	STD02004E
Green	STD02013E
Orange	STD02030E
Red	STD02034E
Violet	STD02039E
White	STD02040E
Yellow	STD02041E

### 1oz

Color:	Item#:
Black	STD02103E
Blue	STD02104E
Green	STD02113E
Orange	STD02130E
Red	STD02134E
Violet	STD02139E
White	STD02140E
Yellow	STD02141E

### 4oz

Color:	Item#:
Black	STD02203E
Blue	STD02204E
Green	STD02213E
Orange	STD02230E
Red	STD02234E
Violet	STD02239E
White	STD02240E
Yellow	STD02241E

### **Mini-Kit: (8x1/4oz) Item#: STD02301P**

*Kit comes with all colors: black, blue, green, orange, red, violet, white, and yellow.*

### **Reg-Kit: (4x1oz) Item#: STD02401P**

*Kit comes with blue, black, red, and green dyes.*

**Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.**

*No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.*