

MasterTech Squeezies

Small, convenient margin marking dyes that can be re-capped for multiple uses.



MasterTech Squeezies are designed to withstand Fixatives, Alcohols, and Clearing Agents. These squeezies are not stains, but special "paints" that coat the tissue area where applied without staining cellular material which could be of diagnostic value.

Identifying a tissue margin is easy as ever with the new MasterTech Squeezies White Marking Dyes. The white dye gives users enhanced visibility during embedding, block orientation, and tissue flotation. Margin appears bright white to the naked eye for fast orientation and dark gray to black under the microscope.

Squeezies come in 1/4 oz. dropper bottle and are available individually, in an eight pack of the same color, or in an assorted color 8 pack, containing one of each.

Individual Squeezies

Color:	Item#:	Color:	Item#:
Black	STD02003E	Red	STD02034E
Blue	STD02004E	Violet	STD02039E
Green	STD02013E	White	STD02040E
Orange	STD02030E	Yellow	STD02041E

Mini-Kit - Same Color/ Pack of 8

Color:	Item#:	Color:	Item#:
Black	STD02003P	Red	STD02034P
Blue	STD02004P	Violet	STD02039P
Green	STD02013P	White	STD02040P
Orange	STD02030P	Yellow	STD02041P



American MasterTech
scientific laboratory supplies



Mini-Kit- All Colors

Black, blue, green, orange, red, violet, white, and yellow.

Item#: STD02301P

Copyright ©2011-2018 American MasterTech Scientific. All rights reserved.

No part of these pages may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, for reasons other than personal use, is strictly prohibited without prior written permission.