

FAST FACTS

MACROSETTE® Processing/Embedding Cassettes

Cassette with lid

Patented



Disposable plastic cassettes designed to hold larger tissue specimens during the embedding process, as well as in a storage cabinet. Dimensions are exactly the same as the ones of a regular HISTOSETTE® cassette but the Macrosette® is twice as high (12mm H.).

Molded from a special high-density polymer, these patented cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage.

The one-piece integral lid eliminates the need for separate steel lids. They can be opened and closed as often as necessary and will always relock securely without danger of specimen loss.

Large labeling areas are located on three sides of the cassettes for your convenience. The anterior writing area is slanted at a 45° angle to make the cassette more suitable to be used with certain types of cassette labeling instruments.

Available in non-cytotoxic, non-metallic color.

Each case contains 3 dispenser boxes of 250 cassettes.

American MasterTech
scientific laboratory supplies

1330 Thurman Street
Lodi, CA 95240

Phone: 800.860.4073
Fax: 209.368.4136

www.americanmastertech.com