

Fast Facts

UNISETTE™ Tissue Processing/ Embedding Cassette

Cassette with lid.



Molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximize fluid exchange and ensure proper drainage. The one-piece integral lid eliminates the need for separate steel lids. The snap-latch and hinge-lock design of the UNISETTE™ prevent early separation of base and lid and allows for one-hand operation. Lids can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Large labeling areas are located on three sides of the cassettes. Anterior writing area is at a 35° angle.

Available in 11 non-cytotoxic, non-metallic colors.

Each case contains 3 dispenser boxes of 500 cassettes.

American MasterTech
scientific laboratory supplies

1330 Thurman Street
Lodi, CA 95240

Phone: 800.860.4073
Fax: 209.368.4136

www.americanmastertech.com